

WREN Laboratories Saliva Collection Kit


For Prescription Use Only
For In Vitro Diagnostic Use Only
For Emergency Use Authorization Only

You have received the WREN Laboratories Saliva Collection Kit. You will self-collect a saliva specimen that will be shipped to WREN Laboratories for COVID-19 testing. Please ship your sample the same day it is collected to ensure timely receipt and processing. Your sample can be shipped any day of the week.

Components of the Saliva Collection Kit:


COLLECTION

- Collection mouthpiece/funnel
- Buffer Cap
- Collection Tube
- Requisition Form
- Collection Instructions



SHIPPING

- Biohazard bag with absorbent pad
- Cardboard box with 2nd absorbent pad
- FedEx envelope
- Pre-paid FedEx return address label
- Shipping Instructions



WARNINGS:

This product has not been FDA cleared or approved;
This product has been authorized by FDA under an EUA for use by the authorized laboratory; This product has been authorized only for the detection of nucleic acid from SARSCoV-2, not for any other viruses or pathogens; and this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the authorization is terminated or revoked sooner.

If you have any questions about your kit, call (203) 208-3464 or email info@wrenlaboratories.com or visit our website at: www.wrencovidtesting.com